

IoT system for Autonomous Realtime Patient Monitoring and Its Application


Unavailable & Espresso Team



Objective

- Using Mobile Sensors
 - Accelerometers
 - Gyroscopes
- Activity classification
 - sleep with turning face up
 - sleep on left's side
 - Sit
 - Fall
- Records
- Alert
 - Idle
 - Fall





Our process

1. Sensors
2. Collect data
3. Cleansing data
4. Model
5. Alert

Sensors




Accelerometer



Gyroscope

Sensors



Collect data

	timestamp	accX	accY	accZ	gyroX	gyroY	gyroZ	pitch	yaw	roll	answer
0	1.637594e+09	0.076492	-0.028275	-1.025864	0.083753	0.020969	-0.069271	0.084725	-0.004381	0.083138	0
1	1.637594e+09	0.069962	-0.116333	-1.000702	0.039438	0.027808	0.023531	0.096635	-0.023273	0.078887	0
2	1.637594e+09	0.066223	-0.074051	-0.969040	-0.024988	0.029181	-0.024695	0.137557	-0.036107	0.099456	0
3	1.637594e+09	0.096603	-0.146362	-0.982025	-0.010441	0.003999	0.007406	0.150668	-0.032262	0.097191	0
4	1.637594e+09	0.100342	-0.147705	-0.982132	-0.012069	-0.001667	0.007258	0.149083	-0.031601	0.101344	0

1. time stamp

2. accX - accelerometers axis-X

3. accY - accelerometers axis-Y

4. accZ - accelerometers axis-Z

5. gyroX - gyroscope axis-X

6. gyroY - gyroscope axis-Y

7. gyroZ - gyroscope axis-Z

8. pitch - gyroscope axis-X

9. yaw - gyroscope axis-Y

10. roll - gyroscope axis-Z

11. answer (y)

Cleansing data

	accX	accY	accZ	pitch	yaw	roll	answer
0	-0.034027	-0.972977	-0.240845	1.393136	0.311162	-0.356315	0
1	-0.032669	-0.995636	-0.138977	1.397252	0.108693	-0.142706	0
2	0.008255	-0.977600	-0.208328	1.362467	-0.072950	0.038964	0
3	-0.080521	-0.966339	-0.223251	1.330597	0.203506	-0.327716	0
4	-0.069138	-0.968109	-0.195251	1.358875	0.180678	-0.302742	0

Explore data

Check mean of each column

	accX	accY	accZ	pitch	yaw	roll
answer						
0	-0.094869	-0.985114	-0.039570	1.421308	1.071202	-1.258685
1	-0.940507	-0.158344	0.290293	0.159052	0.430212	-1.870652
2	-0.226873	-0.168395	-0.944931	0.169737	0.062727	-0.236003



Sit
(class 0)



sleep with turning face up
(class 1)



sleep on left's side
(class 2)

	accX	accY	accZ	pitch	yaw	roll
answer						
0	-0.094869	-0.985114	-0.039570	1.421308	1.071202	-1.258685
1	-0.940507	-0.158344	0.290293	0.159052	0.430212	-1.870652
2	-0.226873	-0.168395	-0.944931	0.169737	0.062727	-0.236003

Model

Splitting Data

```
# Split dataset into training set and test set
X_train, X_test, y_train, y_test = train_test_split(X, y, test_size=0.3, random_state=1) # 70% training and 30% test
```

Building Decision Tree Model

```
# Create Decision Tree classifier object
clf = DecisionTreeClassifier()

# Train Decision Tree Classifier
clf = clf.fit(X_train,y_train)

#Predict the response for test dataset
y_pred = clf.predict(X_test)
```

Evaluating Model

```
print("Accuracy:",metrics.accuracy_score(y_test, y_pred))
```

Accuracy: 1.0

Fall Detection

Fall detection

```
currentSensor = {}  
threshold = 5
```

[+ Code](#)[+ Markdown](#)

```
def fallDetection(accX, accY, accZ, gyroX, gyroY, gyroZ):  
    global currentSensor  
    X = {0:accX, 1:accY, 2:accZ, 3:gyroX, 4:gyroY, 5:gyroZ}  
  
    has_items = bool(currentSensor)  
    if has_items == False:  
        currentSensor = X  
    else:  
  
        diffAccX = abs(accX) - abs(currentSensor[0])  
        diffAccY = abs(accY) - abs(currentSensor[1])  
        diffAccZ = abs(accZ) - abs(currentSensor[2])  
        diffGyroX = abs(gyroX) - abs(currentSensor[3])  
        diffGyroY = abs(gyroY) - abs(currentSensor[4])  
        diffGyroZ = abs(gyroZ) - abs(currentSensor[5])  
  
        print(f'{diffAccX}, {diffAccY}, {diffAccZ}, {diffGyroX}, {diffGyroY}, {diffGyroZ}')  
        if diffAccX > threshold or diffAccY > threshold or diffAccZ > threshold or diffGyroX > threshold or diffGyroY > threshold or diffGyroZ > threshold:  
            alertFallingToMobile()  
  
    currentSensor = X
```



Alert 12

THE TEAM



กรพินธุ์ พิมทาพันธุ์ศรี 6220422056

วัชรวิทย์ ลาพิน 6310422026

วิศัลย์ศยา ทิพย์โรจน์ 6310422027

ธนต์กรณ์ ชื่นบรรลือสุข 6310422031

วริศรา บุญยวง 6310422038

รักษนก ทองคำปิ่น 6310422039

กวิน สิงหเขตต์ 6310422040

ปิยบุญ คุณากรจิตติรักษ์ 6310422047



Thank you!